

rusty-tip

Automating STM/NC-AFM Tip Preparation




Martin Kronberger — github.com/kronberger-droid


Institute for Applied Physics – Surface Physics Group
Technical University Vienna

Introduction	2
Scanning Tunneling Microscope (STM)	3
Non Contact Atomic Force Microscope (ncAFM)	4
Tip Preparation Workflow	5
Software Architecture	7
Outlook	10
References	11

- SPM imaging quality depends critically on the **tip**
- Tip preparation is a manual, experience-driven process
- Goal: develop a software framework in Rust to **automate tip preparation** for STM/ncAFM
- Make the process more **reproducible** and less reliant on individual expertise




- A bias voltage is applied between tip and sample
- Electrons tunnel across the vacuum gap
- The tunneling current I depends exponentially on the tip-sample distance
- Feedback keeps I constant \rightarrow tip traces the surface topography





Non-Contact AFM

- Cantilever vibrates at its resonance frequency
- Tip-surface interaction shifts the frequency by Δf
- Δf is proportional to the force gradient
- A low negative Δf at the surface often indicates a **sharp tip**




Non-Contact AFM

- Cantilever vibrates at its resonance frequency
- Tip-surface interaction shifts the frequency by Δf
- Δf is proportional to the force gradient
- A low negative Δf at the surface often indicates a **sharp tip**



What We Can Automate Now


- **Shape** tip on Au(110), coat with Cu on Cu(110)
- Pulse and reposition until Δf is in the sharp range
- Verify tip stability via bias sweep



Courtesy of D. Kugler

What Comes Next


- Dose O_2 on Cu(110) \rightarrow CuOx rows
- Dip tip into CuOx rows to pick up a cluster
- Image a known surface to verify the tip produces the **expected atomic contrast** ^[1]
- Check tip symmetry by looking at surface impurities



Courtesy of D. Kugler


nanonis-rs ^[3]

- Low-level TCP protocol and parser
- Covers **all** functions available on Nanonis-controlled machines
- High-throughput TCP logger




Action Driver

- Higher-level abstraction over nanonis-rs
- **Action-based system:** each action combines multiple low-level calls
- Makes it possible to log metric and metadata for actions
- Buffered TCP reader for data history analysis



rusty-tip ^[4]

- Uses the Action Driver under the hood
- **State machine** correlating machine states with action loops
- Three states: **Blunt** → **Sharp** → **Stable**
- Logs important metrics for actions



- Automate **tip conditioning** and verification
- Use **image recognition** to evaluate tip quality from images
- Optimize action parameters like pulse voltage via Bayesian optimization
- Extend nanonis-rs to support other SPM workflows beyond tip preparation
- Long-term: fully autonomous experiment — from tip prep to data acquisition

rusty-tip

Automating STM/NC-AFM Tip Preparation



Martin Kronberger — github.com/kronberger-droid

Institute for Applied Physics – Surface Physics Group
Technical University Vienna

1. Mönig et al., “Submolecular imaging by noncontact atomic force microscopy with an oxygen-terminated copper tip,” *Nanoscale* **13**, 18624 (2021).
<https://doi.org/10.1039/d1nr04080d>
2. Hütner et al., “Stoichiometric reconstruction of the $\text{Al}_2\text{O}_3(0001)$ surface,” *Science* **385**, 1241–1244 (2024).
<https://doi.org/10.1126/science.adq4744>
3. nanonis-rs — <https://github.com/kronberger-droid/nanonis-rs>
4. rusty-tip — <https://github.com/kronberger-droid/rusty-tip>